

Flag Tag Agarose

Introduction

TargetMol's Flag Tag agarose gel is prepared by directionally coupling a mouse-derived anti-DYKDDDDK (Flag) monoclonal antibody at high density onto the surface of highly cross-linked agarose beads, enabling specific binding to Flag-tagged proteins.

TargetMol's Flag Tag agarose gel can be used for immunoprecipitation (IP) and co-immunoprecipitation (Co-IP) of proteins and protein complexes, as well as for antibody purification. It is suitable for antigen samples derived from cell lysates, cell culture supernatants, serum, ascites, and other sources.

The product specification refers to the total volume of the suspension.

Product Features

Good physicochemical stability;

Ligand is not easily detached;

High durability;

Easy to use.

Application

1. Immunoprecipitation (IP) and co-immunoprecipitation (Co-IP) of proteins and protein complexes
2. FLAG-tagged protein purification

Product Information

Flag Tag Agarose	Features
Matrix	Highly cross-linked 4% agarose gel
Bead size	30-100 μm
Ligand	Mouse anti-FLAG monoclonal antibody
Binding capacity	≥ 1.0 mg FLAG-tagged protein per mL of gel
Concentration	50% (v/v)

Prepare Reagents

Reagents	Recommended Formulations
Washing Buffer (1 \times)	TBST: 50 mM Tris-HCl, 150 mM NaCl, 0.1%(v/v) Tween-20, pH7.4
Flag Peptide Elution Buffer	PBS, 1 mg/mL 3 \times Flag peptide (TP1274), pH 7.4
Acidity Elution Buffer	0.1 M Glycine, 0.1% (v/v) Tween-20, pH2.5
Neutralization Buffer	1 M Tris-HCl, pH 9.0

Instructions

1. Preparation of Cell Lysate

Select an appropriate lysis buffer to treat the cell samples. Prepare the cell lysate according to standard procedures. Keep on ice for immediate use, or store at -20 °C for long-term preservation.

2. Gel Pre-treatment

Vortex the agarose gel for 1 min to fully resuspend it. Transfer 25-50 μL of the gel suspension into a 1.5 mL EP tube. Add 500 μL of Washing Buffer to the tube and gently invert several times to resuspend the gel.

Centrifuge the tube at 1,000 rpm for 5 min. After the gel has settled at the bottom, discard the supernatant. Repeat the washing step once.

3. Immunoprecipitation Procedure

Add the prepared sample containing FLAG-tagged protein into the EP tube. Place the tube on a rotator and incubate at room temperature for 1-2 h, or at 4 °C for 2-4 h.

After incubation, centrifuge at 1,000 rpm for 5 min. Collect the gel and discard or retain the supernatant for further analysis.

Add 1,000 µL of Washing Buffer to the tube, and gently mix for 5-10 min. Centrifuge at 1,000 rpm for 5 min, collect the gel, and discard the supernatant. Repeat the wash step two more times.

4. Antigen Elution

The following antigen elution methods are provided. Users can choose the appropriate method based on downstream applications.

A. Denaturing Elution Method:

The eluted sample is suitable for SDS-PAGE analysis. Add 100 µL of SDS-PAGE Loading Buffer (self-prepared) to the EP tube and mix thoroughly. Heat at 95 °C for 5 min. Then centrifuge to separate the gel, collect the supernatant, and proceed with SDS-PAGE analysis.

B. Neutral Elution Method:

Add 50 µL of Flag Peptide Elution Buffer to the EP tube. Incubate on a rotator at 37 °C for 5-10 min (extend incubation time if below 37 °C). Then centrifuge to separate the gel and collect the supernatant. To improve antigen recovery, the elution step can be repeated.

C. Acidic Elution Method:

Add 100 µL of Acidity Elution Buffer to the EP tube. Incubate on a rotator at 37 °C for 5-10 min. Then centrifuge to separate the gel and collect the supernatant. If neutralization is required, add 50 µL of Neutralization Buffer to 100 µL of the eluate to adjust the pH to neutral.

Storage

Store at 4 °C for 2 years.

Precautions

1. The gel should be stored in the storage solution to prevent drying.
2. Before taking the agarose gel from the storage tube, mix thoroughly to ensure a uniform suspension. Avoid introducing bubbles during handling.
3. Depending on experimental needs, the operator may use the supernatants collected during the antibody-binding and antigen-binding steps to assess the binding efficiency between antibody, antigen, and gel.
4. In IP experiments, the binding affinity between different antibodies and antigens may vary. If the buffer system provided in this kit does not yield satisfactory results, the operator may optimize or prepare customized buffers as needed.
5. This product is for R&D use only, not for diagnostic procedures, food, drug, household, or other uses.
6. It's advisable to wear a lab coat and disposable glove.

TargetMol US

www.targetmol.com sales@targetmol.com 1-781-999-5354

36 Washington Street, Wellesley Hills, MA 02481 USA

TargetMol EU

www.targetmol.com sales@targetmol.com +43(0)676/7860258

Hafenstraße 47-51, 4020 Linz, Austria



LinkedIn



Facebook



TargetMol Email