

Mouse CD3/CD28 T Cell Activation Beads

Introduction

TargetMol's Mouse CD3/CD28 T Cell Activation Beads are a convenient and efficient tool for activating and expanding T cells without the need for feeder cells (such as antigen-presenting cells) or antigens. The beads are 5 μm in diameter, made of inert superparamagnetic material, and are comparable in size to antigen-presenting cells. They are covalently conjugated with anti-CD3 and anti-CD28 antibodies. In cell culture, the addition of recombinant mouse IL-2 can further promote the expansion of the T cell population. After activation and expansion, the beads can be easily removed using a magnetic separator.

Related Products

1. Mouse Cells

	Spleen	Lymph Nodes	Peripheral Blood	Bone marrow	Tumor tissue
CD3 ⁺ T Cell	C0061		/	/	/
CD4 ⁺ Cell	C0062 (Preferred) , C0067 (Optional)		C0067	/	C0067
CD8 ⁺ Cell	C0063 (Preferred) , C0068 (Optional)		C0068	/	C0068
Neutrophils	C0064	/	C0064	C0064	/
CD3 ⁺ Cell Depletion	C0152	C0152	/	/	/
CD3/CD28 T Cell Activation	C0180	/	/	/	/

2. Human Cells

	Peripheral Blood	Cord Blood
CD3 ⁺ T Cell	C0065	/
CD34 ⁺ Cell Enrichment	C0066	C0066
CD4 ⁺ T Cell	C0148	/
CD8 ⁺ T Cell	C0149	/
CD3/CD28 T Cell Activation	C0150	/
CD66b ⁺ Cell	C0151	/

Features

1. High Efficiency: Rapidly activates and expands T cells without the need for antigen-presenting cells or antigens.
2. High Viability: Beads are non-cytotoxic, ensuring T cells maintain strong vitality and functionality.
3. Easy to Use: Enables target cell separation with a magnetic stand.

Application

1. Intended for activation of mouse T cells, e.g. CD3⁺ T cells, CD4⁺ T cells, or CD8⁺ T cells.
2. The activated T cells can be analyzed shortly after activation (for trans-fec-tion/transduction or to study T cell receptor signaling, proteomics, or gene expression). T cells can be left in culture to differentiate into T helper cell subsets, T cell proliferation or expansion of polyclonal T cells.

Product Information

Mouse CD3/CD28 T Cell Activation Beads	Features
Concentration	1 \times 10 ⁸ beads/mL
Solution	PBS pH 7.4 (0.1% BSA, 0.1% proclin-300)

Prepare Reagents

1. **Washing buffer:** Use PBS containing 2 mM EDTA and 2% fetal bovine serum (FBS), or PBS containing 2 mM EDTA and 0.5% BSA. Sterilize by passing through a 0.22 μ m filter before use.
2. **Cell culture medium:** Use RPMI 1640 supplemented with 2 mM L-glutamine and 10% FBS (optionally with 100 U/mL penicillin-streptomycin). For T cell expansion, recombinant mouse IL-2 (e.g., TMPY-02788) can be added at 30 U/mL. Alternatively, commercially available T cell culture media can be used, with growth factors in the formulation adjusted flexibly according to experimental needs.

Instructions

Take isolated mouse CD3+ T Cells:

Wash Magnetic Beads

1. Resuspend the magnetic beads in the vial (i.e. vortex for >30 sec, or tilt and rotate for 5 min).
2. Transfer the desired volume of magnetic beads to a tube. Add an equal volume of Buffer (If the bead volume is less than 1 mL, add wash buffer to bring the total volume to 1 mL). Gently vortex for 5 sec or pipette up and down to mix. Avoid generating bubbles during handling.

Note: cell culture medium can be used instead of wash buffer.

3. Place the tube on a magnet for 3 min and discard the supernatant.
4. Wash the beads again following steps 2 and 3, this time using cell culture medium as the wash solution. In total, wash the beads twice.
5. Resuspend the washed beads in cell culture medium. For example, if 25 μ L of beads are washed, finally resuspend them in 1 mL of cell culture medium.

Activate Mouse T Cells (Take a 96-well plate as an example)

1. Adjust the T cell concentration to 1×10^7 cells/mL, add 25 μ L of the cell suspension to each well (2.5×10^5 T cells per well), and supplement with culture medium to a final volume of 100 μ L.
2. Add 100 μ L pre-washed and resuspended beads to obtain a bead-to-cell ratio of 1:1.
3. Incubate in a humidified 5% CO₂ incubator at 37°C, according to your specific experimental requirements.
4. Harvest the activated T cells collected in 24-48 h and use directly for further analysis.
5. For flow cytometry applications, remove the beads prior to staining. Place the tube on a magnet for 1 – 2 min to separate the beads from the solution. Transfer the supernatant containing the cells to a new tube.

Note: To increase recovery of T cells in the supernatant, resuspend the bead/cell suspension thoroughly by pipetting prior to magnetic separation (some cells bind strongly to the beads upon activation and for 2 – 3 days thereafter). To further increase cell recovery, culture the bead-bound cell fraction overnight and repeat step 5.

Expand Mouse T Cells

1. Check the activation status of T cell after 48 hour. If the cells are in good condition, gently pipette to resuspend them and proceed with passaging. For cells with a slower proliferation rate, perform a half-medium change first: carefully remove 100 μ L of supernatant from each well, add an equal volume of medium, and gently pipette to mix. Continue culturing for another 1-2 days before passaging. Finally, add 30 U/mL recombinant mouse IL-2 to the culture medium to support subsequent T cell proliferation and activation.
2. Incubate in a humidified 5% CO₂ incubator at 37 °C, according to your specific experimental requirements.
3. Examine cultures daily, noting cell size and shape. Cell shrinking and reduced proliferation rate is typically observed in exhausted cell cultures.
4. It is recommended not to handle the cells during the initial activation period (first 2 days). After 2 days, observe them regularly. When the culture medium turns yellow or the cell density becomes too high, perform medium change or passaging. Supplement with 30 U/mL recombinant mouse IL-2 each time.
5. Count the cells periodically. When the cell density exceeds 2.5×10^6 cells/mL, split cultures back to a density of $0.5-1 \times 10^6$ cells/mL.

Storage

Store at 4 °C for 2 years.

Precautions

1. Avoid freezing the components of the kit. Store in the storage buffer to prevent drying.
2. Before withdrawing the beads, shake thoroughly to ensure a uniform suspension. Avoid generating bubbles during handling.
3. Use high-quality tips and test tubes to avoid sample loss due to adhesion.
4. The product is for R&D use only, not for diagnostic procedures, food, drug, household or other uses.
5. Please wear a lab coat and disposable gloves.

TargetMol US

 www.targetmol.com  sales@targetmol.com  1-781-999-5354
 36 Washington Street, Wellesley Hills, MA 02481 USA

TargetMol EU

 www.targetmol.com  sales@targetmol.com  +43(0)676/7860258
 Hafenstraße 47-51, 4020 Linz, Austria



LinkedIn



Facebook



TargetMol Email