

### 1.25% Avertin (Anesthesia For Mice)

---

#### Description

1.25% Avertin (Anesthesia For Mice) is a sterile solution optimized for short-acting anesthesia in laboratory animals, particularly mice. It uses high-purity tribromoethanol as the main active ingredient, with a small amount of tert-amyl alcohol as a solubilizing agent, and is formulated with physiological saline to a standard concentration of 1.25%.

Tribromoethanol is a central nervous system depressant that enhances GABAA receptor-mediated chloride ion influx, suppressing neuronal depolarization and thereby inhibiting neural excitation. This rapidly induces sedation and anesthesia. As a fast-acting anesthetic with a short duration, Avertin is widely used in acute mouse experiments, such as surgeries, cardiac puncture, perfusion, and sample collection. The 1.25% concentration is the commonly recommended dose for mice, balancing anesthesia depth and safety.

This product is ready to use that has been sterilized by filtration to ensure safety and stability. It is suitable for direct intraperitoneal injection for laboratory animal anesthesia, with no need to dissolve or filter.

#### Features

1. **Ready-to-use:** use directly and avoid exposure to toxic gases.
2. **High-purity raw materials:** Made with analytical-grade tribromoethanol, sterile saline, and solubilizing agents.
3. **No equipment required:** No need for a vaporizer or anesthesia machine, suitable for standard laboratory environments.
4. **Strict sterility:** Prepared by 0.22 µm sterile filtration, suitable for perfusion and other high-cleanliness experiments.
5. **Stable packaging:** Light-protected, cold-resistant storage bottles to ensure product activity and safe use.

#### Application

For use in rodent (mice) anesthesia only.

#### Instructions

**Administration route:** Intraperitoneal injection (i.p. injection)

**Recommended dosage:** Approximately 200-300 µL/10 g mouse body weight

**Onset of anesthesia:** within 5 minutes

**Duration of anesthesia:** 20–90 minutes (varies depending on dose and animal condition)

Note: Please adjust the dosage based on the animal's condition and experimental requirements to avoid overdose.

#### Storage

Store at 4°C protected from light; valid for one year.

#### Precautions

1. Different strains and ages of animals require different doses of anesthetics. It is recommended to start with a low dose and conduct preliminary experiments to optimize the final dosing.
  2. Use as soon as possible after opening. Prolonged storage may cause precipitation or decomposition of components, resulting in reduced anesthetic efficacy and even peritonitis. Allow the solution to return to room temperature before use. Cloudiness or crystal formation is normal and can be dissolved in a 37°C water bath.
  3. Within the recommended range, the larger the dose of Avertin, the longer the duration of deep anesthesia and recovery time. A heating pad or heat lamp can help shorten recovery time.
-

4. Avertin is suitable for short-term survival or non-survival surgeries and terminal euthanasia procedures. Adjust the initial dose according to the duration of the procedure, and if necessary, administer approximately one-third of the original anesthetic dose to extend the anesthesia duration.
5. Suitable for use in pregnant mice; does not affect pregnancy and has minimal impact on fetuses.
6. Do not use Avertin repeatedly on the same animal, as it may cause peritonitis and death. If the animal has been anesthetized with Avertin and has recovered, alternative methods (e.g., isoflurane) should be used for subsequent anesthesia.
7. The product is for R&D use only, not for diagnostic procedures, food, drug, household or other uses.
8. Please wear a lab coat and disposable gloves.

## TargetMol US

 [www.targetmol.com](http://www.targetmol.com)  [sales@targetmol.com](mailto:sales@targetmol.com)  1-781-999-5354  
 36 Washington Street, Wellesley Hills, MA 02481 USA

## TargetMol EU

 [www.targetmol.com](http://www.targetmol.com)  [sales@targetmol.com](mailto:sales@targetmol.com)  +43(0)676/7860258  
 Hafenstraße 47-51, 4020 Linz, Austria



LinkedIn



Facebook



TargetMol Email