

0.25% Trypsin-EDTA, phenol red (1×)

Description

TargetMol's 0.25% Trypsin-EDTA, Phenol Red (1×) is a ready-to-use cell dissociation solution specifically designed for routine passaging, harvesting, and single-cell preparation of adherent cells. The formulation contains trypsin and EDTA, which effectively disrupt cell-substrate and cell-cell interactions, enabling rapid and efficient cell detachment. The inclusion of phenol red serves as a pH indicator, allowing real-time monitoring of the dissociation process. This optimized formulation ensures efficient cell detachment while maximally preserving cell membrane integrity and biological activity, making it suitable for a wide range of mammalian adherent cell types. Supplied as a 1× working solution, the product can be used directly without dilution. It is ideal for routine cell culture applications, cell passaging, and pre-freezing preparation, and is an essential, commonly used reagent in cell culture laboratories.

Features

- **Ready to Use:** No additional dilution or preparation required; can be used directly for cell dissociation, saving experimental setup time.
- **Strong & Fast Dissociation:** 0.25% trypsin efficiently dissociates tightly adherent cells, significantly reducing dissociation time.
- **Synergistic Action with EDTA:** EDTA chelates Ca^{2+} and Mg^{2+} , enhancing the disruption of cell-cell junctions and improving dissociation uniformity.
- **Phenol Red Indicator Included:** Allows real-time monitoring of the dissociation process, with intuitive color changes indicating progress.
- **Consistent Quality:** Manufactured with high-purity trypsin, sterile-filtered to ensure cell culture safety and reproducibility.
- **Broad Compatibility:** Suitable for routine passaging and dissociation of a wide range of mammalian cell lines.

Application

For the dissociation, passaging, and collection of adherent cells.

Instructions

1. Preparation

a) Take the trypsin solution from the refrigerator and prewarm it at 37 °C in a water bath or incubator for 5-10 minutes.

Note: Do not prewarm the entire bottle; only take out the volume you need.

b) Remove the old medium from the culture dish or flask, and gently rinse once with sterile PBS, Hanks, or serum-free medium to remove residual serum.

Note: Serum inhibits trypsin activity.

2. Trypsinization

a) Add an appropriate amount of trypsin solution to just cover the cell layer. Incubate at room temperature for 1-5 minutes, gently swirling the culture vessel occasionally to ensure even contact with the cells.

b) Observe under a microscope. When the cells begin to round up and detach from the bottom, dissociation is sufficient.

3. Termination of Trypsinization

a) Immediately add an equal or double volume of complete medium containing serum to stop the trypsin activity.

b) Gently pipette to mix, ensuring the cells are fully detached into a single-cell suspension.

4. Cell Collection and Culture

a) Transfer the cell suspension to a centrifuge tube and centrifuge at 1,000 rpm for 3-5 minutes.

b) Discard the supernatant and resuspend the cells in fresh medium. Proceed with counting, passaging, or seeding as required by your experiment.

Storage

Store at 4 °C for 3 months; -20 °C for 2 years.

Precautions

1. After use, immediately seal and store at 4 °C. For long-term storage, keep at -20 °C and avoid repeated freeze-thaw cycles.

2. Avoid excessively long trypsin dissociation, as this may cause cell damage or death.
3. For sensitive cells (e.g., primary cells or stem cells), consider shortening dissociation time or diluting the trypsin concentration.
4. This product contains phenol red as a pH indicator. If the pH slightly decreases for any reason, the solution may turn from red to orange. The orange solution can still be used normally, or, if needed, adjust the pH slightly with sterile 2 M NaOH before use.
5. Observe aseptic techniques to prevent microbial contamination.
6. The product is for R&D use only, not for diagnostic procedures, food, drug, household or other uses.
7. This product may irritate skin, eyes, and the respiratory tract. Please wear a lab coat and disposable gloves.

How to Select Trypsin Cell Dissociation Solutions

	Cells are sensitive to trypsin & dissociation time is difficult to control	Strong dissociation	To monitor the dissociation process	Without Phenol Red	Without EDTA
C0200 0.05% Trypsin-EDTA, phenol red (1x)	✓		✓		
C0201 0.25% Trypsin-EDTA, phenol red (1x)		✓	✓		
C0202 0.25% Trypsin-EDTA (1x)		✓		✓	
C0203 0.25% Trypsin, phenol red (1x)			✓		✓

TargetMol US

 www.targetmol.com
 sales@targetmol.com
 1-781-999-5354
 36 Washington Street, Wellesley Hills, MA 02481 USA

TargetMol EU

 www.targetmol.com
 sales@targetmol.com
 +43(0)676/7860258
 Hafenstraße 47-51, 4020 Linz, Austria



LinkedIn



Facebook



TargetMol Email